

TELECINEMA AND SCAN BY EDL

Telecine
from the negative
while preparing
the reels
to be able to pull the
used images in the end



- Punching every lab roll



6 - 10 € / min



- Creating a file for each lab roll starting with the punch.

- Giving time-code 01:00:00:00 to the punch of the first lab roll, and the following hour for each of the remaining ones.

EDITING
DIGITALLY



- Exporting an EDL that allow scanning only the used images in the lab roll.



SCAN BY EDL

50 - 200 € / min



Color timing
in digital